

Solubility of Riboflavin for Coverage Testing

by Chris Bautz

CIP Summit
September 18, 2019



© 2019 by Electrol Specialties Company

Proud to be an employee owned company

Riboflavin Spray Coverage Testing

- Establishes that all vessel process surfaces are contacted
- Riboflavin fluoresces under UV light for visual confirmation
- No visible fluorescence after test required for passing
- Not a cleaning test – evaluates coverage only!
- Sounds simple, but it's a bit more complicated than that

Riboflavin Testing – False Failures

- A failed coverage test can be caused by factors other than flaws in the design of the spray devices or patterns
 - Trapping riboflavin in manway gasket or other seal and releasing back into vessel after opening
 - Contaminants or soils on the surface may prevent proper rinsing
- Incompletely dissolved riboflavin can give false results also!
- Unpredictability of results can be frustrating

Clouding the Solution



Riboflavin Mixture vs Riboflavin Solution

Proud to be an employee owned company

Riboflavin Testing – False Failures

- Undissolved riboflavin powder can stick to interior tank surfaces and prevent proper rinsing
- Droplets “reconcentrate” as remaining riboflavin dissolves
- Riboflavin on the spray devices themselves?!?
 - Not unheard of
 - High likelihood of being a false failure
 - Especially with straight spray devices

False Failure Illustrated on Sprayball



Initial inspection after 3 x 60 sec
rinse bursts: Pass!




5 Minutes Later: Wait, it failed!?

Proud to be an employee owned company

Riboflavin Solution Preparation


- Preparation is Important
 - BPE specifies a riboflavin solution, not mixture
 - Riboflavin is harder to dissolve in water than you may think
 - “Slightly soluble in water” according to riboflavin manufacturer literature
- Heating the water helps



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps

- 0.2 g/L riboflavin aqueous solution
- Solubility tested at range of temps
 - $\approx 25^{\circ}\text{C}$ (77°F)
 - $\approx 40^{\circ}\text{C}$ (104°F)
 - $\approx 50^{\circ}\text{C}$ (122°F)
 - $\approx 60^{\circ}\text{C}$ (140°F)
 - $\approx 70^{\circ}\text{C}$ (158°F)



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps

- Riboflavin in RO water at 25° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Obvious cloudiness due to incompletely dissolved granules



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps


- Riboflavin in RO water at 40° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Less cloudy, but ribo obviously still not completely dissolved



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps


- Riboflavin in RO water at 50° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Less cloudy still
- Stir bar becoming visible, but undissolved powder still present



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps


- Riboflavin in RO water at 60° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Small amount of riboflavin powder still visible when mixing stopped



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps


- Riboflavin in RO water at 70° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Riboflavin particles finally appear to be completely dissolved



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps


- RO water on 16AMP solid end cap
- <20 Ra mechanical
- 10x magnification



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps


- Riboflavin in RO water at 25° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Obvious cloudiness due to riboflavin powder being incompletely dissolved



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps


- Riboflavin in RO water at 40° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Less cloudy, but ribo obviously not completely dissolved



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps


- Riboflavin in RO water at 50° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Less cloudy still
- Stir bar becoming visible, but undissolved powder still present



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps


- Riboflavin in RO water at 60° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Small amount of riboflavin powder still visible when mixing stopped



Proud to be an employee owned company

Dissolving Riboflavin at Different Temps

- Riboflavin in RO water at 70° C
- Mixed for 5 minutes on high setting of magnetic lab mixer
- Riboflavin particles finally appear to be completely dissolved



Proud to be an employee owned company

ESC Riboflavin Mixing Recommendations

- Heat at least 50% of final solution volume of water to 80°C
- Mix the riboflavin powder into the hot water with vigorous agitation until no undissolved powder remains
- Add the remaining water to the hot riboflavin solution to achieve final solution volume & concentration
- Mix vigorously and allow to cool before use

Proud to be an employee owned company